

Product Application

contact techserv@promega.com

Isolating DNA from Cotton Leaf Tissue Using the ReliaPrep[™] gDNA Tissue Miniprep System

Isolate high quality, amplifiable DNA from cotton leaf tissue using the ReliaPrep[™] gDNA Tissue Miniprep System.

Ki+·	PoliaPron™ aDNA Tissue Minipron System (Cat. #	
Kit.	A2051)	This protocol was developed by Promega Applications Scientists and is intended
Sample Type:	Cotton leaf tissue	for research use only. Users are responsible for determining
Input:	Up to 25mg	suitability of the protocol for their application.
Materials Required: •	ReliaPrep™ gDNA Tissue Miniprep System (Cat. # A2051) 2.0ml ceneus ten tubes	Further information can be found in Technical Manual #TM345, available at: <u>www.promega.com/protocols</u> OR For further information, please

- 2.0mL screw-top tubes
- Homogenization steel bead
- Bead-beating device (e.g., MP Biomedicals FastPrep[®]-24 Instrument)
- Microcentrifuge

Protocol:

- 1. Using a 5mm punch, place desired number of punches (up to 25mg) into a 2ml screw-top tube.
- 2. To each sample add:
 - 100µl of Tail Lysis Buffer (TLA)
 - 300µl of Cell Lysis Buffer (CLD)
 - 20µl of RNase A Solution
 - 20µl of Proteinase K
- 3. Using the bead-beating device, homogenize samples for desired time (e.g., FastPrep[®]-24 Instrument at 6m/s, 40 seconds).
- 4. Centrifuge samples in a microcentrifuge at max speed for 1 minute.
- 5. Incubate at room temperature for 10 minutes.
- 6. Centrifuge samples at max speed for 1 minute to reduce foaming.
- 7. Add 250µl of Binding Buffer (BBA) to each sample and vortex for 10 seconds.
- 8. Centrifuge samples at max speed for 2 minutes and transfer liquid supernatant to a ReliaPrep[™] Binding Column inside a collection tube.
- 9. Centrifuge samples at max speed for 1 minute. Transfer column to a new collection tube, discard the flow through and used collection tube.
- 10. Add 500µl of Column Wash Solution (CWD) to the sample and centrifuge at max speed for 2 minutes. Repeat this wash step for a total of 3 times, discarding liquid and collection tubes after every wash.
- 11. Place the column in a clean 1.5ml microcentrifuge tube.



12. Add 50µl Nuclease-Free Water to the column. Centrifuge for 1 minute at max speed. Eluates are ready for use in downstream applications.

Results:

Table 1. Cotton leaf DNA concentrations, yields, and purity based on quantitation using theQuantiFluor® ONE dsDNA System (Cat. # E4871) and the NanoDrop-1000. DNA of high purity wasrecovered with purity ratios for samples >1.7. N=3.





Figure 1. Inhibition analysis of purified cotton leaf DNA. DNA samples were serial diluted 1:10 and 10:100. For a sample diluted 10-fold, Δ Cq values are expected to be 3.3. Δ Cq values significantly less than 3.3 may indicate the presence of inhibitors. Δ Cq values of DNA from cotton leaf samples indicate little to no inhibition of the serial diluted eluates. N=3.