

or by e-mailing Technical Services techserv@promega.com

Clean-Up of In Vitro Transcription Reactions Using the ReliaPrep™ **RNA Clean-Up and Concentration System**

Clean-up in vitro transcription reaction rapidly and without phenol chloroform extraction

Kit. F	ReliaPren™ RNA Clean-Un and Concentration System (71073)	
Analyses:	NanoDrop [™] , QuantiFluor RNA System, TapeStation	This protocol was developed by Promega Applications Scientists and is intended for research use only.
Sample Type(s):	In vitro transcription reactions	Users are responsible for determining suitability of the protocol for their
Input:	20µL	application.
Materials Requir	ed:	Technical Manual #TB298, available at: www.promega.com/protocols

- Microcentrifuge
- Vortex
- 100% Isopropanol
- 95-100% Ethanol

Protocol:

- 1. Pipet 20µL of in *in vitro* transcription reactions into a 1.5mL microcentrifuge tube.
- 2. Add 12.5µL of Membrane Binding Solution and vortex for 5 seconds.
- 3. Add 37.5µL of 100% isopropanol.
- 4. Load RNA sample onto a ReliaPrep[™] Minicolumn seated in a Collection Tube and centrifuge for 30 seconds.
- 5. Remove column, and discard the contents of the Collection Tube. Reseat the column into the same Collection Tube.
- 6. Add 200µL of Column Wash Solution (CWE) and centrifuge for 15 seconds. Remove column, and discard the contents of the Collection Tube. Reseat the column into the same Collection Tube.
- 7. Wash with 300µL of RNA Wash Solution (RWA) and centrifuge for 15 seconds. Repeat wash with 300µL of RNA Wash Solution (RWA) and centrifuge again.
- 8. Remove column and discard the contents of the collection tube. Reseat the column into the same collection tube, centrifuge for 1 minute to dry the column and then transfer the column to an Elution Tube.
- 9. Pipet 15µL of Nuclease-Free Water to the center of the column matrix and centrifuge for 30 seconds.
- 10. For maximum recovery, repeat elution with an additional 15µL of Nuclease-Free Water.



Product Application

Results:

RNA was transcribed from the pGEM[®] Express Positive Control Template using the T7 RiboMax[™] Express Large-Scale RNA Production System following the standard protocol (TB298). Reactions were purified using phenol chloroform extraction as recommended by the technical manual or the ReliaPrep[™] RNA Clean-Up and Concentration System (TM541). One of the phenol chloroform extractions was not included in this data set due to phenol contamination.



Figure 1. DNA yield measured by absorbance using NanoDrop[™] (left) and QuantiFluor[®] RNA System (right) for *in vitro* transcription reactions purified using phenol chloroform extraction or the **ReliaPrep[™] Clean-Up and Concentration System.** Mean ± Standard Deviation of n=3 for phenol chloroform extraction and n=4 for ReliaPrep[™] RNA Clean-Up and Concentration System are shown.



Figure 2. RNA integrity for *in vitro* transcription reactions purified using phenol chloroform extraction (A-C) or the ReliaPrep[™] Clean-Up and Concentration System (D-F). RNA integrity was tested using RNA ScreenTape on an Agilent 4200 TapeStation according to the technical manual as shown above.