

Product Application

Coffee (Coffea arabica) Leaf DNA Purification

Isolate high quality, amplifiable DNA from coffee leaf tissue using the Maxwell[®] 16 System.

Kit:	Maxwell [®] 16 LEV Plant DNA Kit (Cat. #AS1420)	
Analyses:	GoTaq [®] qPCR, QuantiFluor [®] quantitation, gel	This protocol was developed by Promega Applications Scientists and is intended for research use only.
Sample Type(s):	Coffea arabica young plant leaves	Users are responsible for determining suitability of the protocol for their application.
Input:	1-5 punches of leaf tissue, 1 punch = ~ 2.5mg	
Materials Required:	 Maxwell[®] 16 Instrument (Cat. #AS2000) with firmware version 4.97 or later Maxwell[®] 16 LEV Plant DNA Kit (Cat. #AS1420) 2.0mL screw-top tubes Homogenization steel bead 	Further information can be found in Technical Manual #TM414, available at: www.promega.com/protocols

- Bead-beating device (e.g., MP Biomedicals FastPrep®-24 Instrument)
- Microcentrifuge

Protocol:

- 1. Using a 5mm punch, place desired number of punches (up to 5 punches) into a 2ml screw-top tube.
- 2. Add 300µl of Tissue Lysis Buffer (TLA) to each sample tube.
- 3. Using the bead-beating device, homogenize samples for desired time (e.g., FastPrep®-24 Instrument at 4M/S, 20seconds x 4 times with 20 second delay between each time).
- 4. Centrifuge samples in a microcentrifuge at max speed for 2 minutes.
- 5. Add 300µl of Nuclease Free Water to well #1 of each Maxwell cartridge.
- 6. Transfer the entire volume of supernatant to well #1 of the Maxwell cartridge.
- 7. Place one of the supplied elution tubes into the sample rack and add 50µl of the supplied Elution Buffer for each sample.
- 8. Place the plunger in the indicated position of the cartridge.
- 9. Select LEV configuration on the Maxwell[®] and select method: RUN, DNA: Plant. Start run.



Product Application

Results

The above protocol was tested using 1-5, 5mm punches of Arabica coffee plant leaves.





M=Benchtop 1kb ladder

Coffee leaf DNA yield and quality. **Left Panel:** DNA yield from 1-5 punches was quantified using the QuantiFluor[®] ONE dsDNA System (Cat. # E4871). **Right Panel:** Examples of high molecular weight DNA purified from the indicated number of punches and analyzed on a 1.0% Agarose gel with 5µl DNA eluate per lane. M = BenchTop 1kb DNA Ladder (Cat. # G7541).



Performance of purified DNA in qPCR. Left Panel: Analysis of purified DNA using GoTaq[®] qPCR Master Mix (Cat. # A6001) and Plant Universal Primers (1), using 1µl DNA eluate per 25µl reaction. **Right Panel:** Δ Ct of serially diluted samples indicates minimal qPCR inhibition.

References:

1. Wang, J., *et al.*, (2011). Universal endogenous gene controls for bisulfite conversion in analysis of plant DNA methylation. *Plant Methods*. **7**, 39.