

# **Product Application**

## Rice (Oryza sativa) Kernel DNA Purification

Isolate high quality, amplifiable DNA from brown or white rice kernels using the Maxwell<sup>®</sup> 16 System.

Kit:	Maxwell <sup>®</sup> 16 FFS Nucleic Acid Extraction Kit	
	(Cat. #X9431)	This protocol was developed by Promega Applications Scientists and is intended
Analyses:	GoTaq <sup>®</sup> qPCR, QuantiFluor <sup>®</sup> quantitation, gel	for research use only.
Sample Type(s):	<i>Oryza sativa</i> kernels (brown or white medium grain or basmati rice)	Users are responsible for determining suitability of the protocol for their application.
Input:	up to 200mg ground or whole kernels	For further information, please contact techserv@promega.com
Materials Required:		

- Maxwell<sup>®</sup> 16 Instrument (Cat.# AS2000) with firmware version 4.97 or later
- Maxwell<sup>®</sup> 16 FFS Nucleic Acid Extraction Kit (custom Cat.# X9431)
- <u>Optional</u>: RNase A Solution (Cat.# A7973)
- 2.0mL screw-top tubes
- CTAB Buffer: 2% CTAB, 1.4M NaCl, 0.1M Tris 10mM EDTA pH 8.0
- <u>Homogenization protocol only</u>: Coffee grinder or other grinding apparatus
- Microcentrifuge
- Heat block 65°C

#### Protocol:

<u>Ground kernels</u> (Recommended for brown rice): Grind kernels to a fine powder in liquid nitrogen with a mortar and pestle or in a coffee grinder with pulsing. Weigh up to 200mg of ground kernel into a 2ml tube.

Whole kernels: Weigh up to 200mg of whole kernels (about 10 kernels) into a 2ml tube.

- 1. Add 1ml CTAB buffer spiked with 40μl Proteinase K and <u>optionally</u> 20μl RNase A Solution to prepared rice kernels.
- 2. Vortex vigorously and incubate 90 minutes at 65°C.
- 3. Centrifuge 10min at 16,000 x g.
- 4. Place the cartridge into the Maxwell<sup>®</sup> LEV cartridge rack and remove the seal.
- 5. Pipet 300µl sample supernatant and 300µl Lysis Buffer into well #1 of the Maxwell<sup>®</sup> 16 cartridge.
- 6. Place one of the supplied elution tubes into the sample rack and add 100µl of the supplied elution buffer.
- 7. Select LEV configuration of the Maxwell<sup>®</sup> Instrument and select method as follows: RUN, DNA: Plant. Start Run.



# **Product Application**

#### Results

Medium Grain Rice

Rice Type	Input	[DNA] (ng/μl)
White	Kernel	35 ±10.4
	Ground	17.5 ±2.1
Brown	Kernel	2.1 ±0.7
	Ground	13 ±2.8

DNA concentration and yield. Left Panel:

Concentration of DNA extracted from 100mg white and brown medium grain rice whole or ground kernels and eluted in 100µl.

Concentration was calculated using the



QuantiFluor<sup>®</sup> ONE dsDNA System (Cat. # E4871). **Right Panel:** Yield of amplifiable DNA determined by real-time PCR using GoTaq<sup>®</sup> qPCR Master Mix (Cat.# A6001), rice DNA specific primers, and a Rice Genomic DNA standard (Zyagen, Cat. #PLG-1004).

### **Basmati Rice**

Rice Type	Input	[DNA] (ng/µl)
White whole kernel	100mg	28.3 ±5.1
Brown ground	100mg	23.7 ±11.5

**DNA concentration and yield: Left Panel**: Concentration of DNA extracted from white kernels or brown ground Basmati rice and eluted in 100µl. Concentration was calculated using the QuantiFluor<sup>®</sup> ONE dsDNA System (Cat.# # E4871). **Right Panel:** Yield of amplifiable DNA determined by real-time PCR using GoTaq<sup>®</sup> qPCR Master Mix (Cat.# A6001), rice DNA



specific primers, and a Rice Genomic DNA standard (Zyagen, Cat. #PLG-1004).

### Single Rice Kernel

To isolate DNA from single rice kernels, minor protocol modifications are suggested to maximize DNA concentration: Use half the CTAB and proteinase K volumes. After heating, add the entire sample (~520µl) to well #1, and elute in 50µl.



PCR amplification of extracted rice kernel DNA. The gel image

shows PCR products from DNA extracted from a single rice kernel and amplified using the Rice Multiple-X<sup>™</sup> PCR Kit (SolGent Co., Ltd); 5µl gDNA was added to a 25µl reaction. Six different types of rice were tested for SNP identification.