

Product Application

ccfDNA extraction from plasma, urine and cerebrospinal fluid using Maxwell[®] RSC System

Isolate amplifiable ccfDNA from plasma, urine, and cerebrospinal fluid using the Maxwell® RSC System

Kit·	Maxwell [®] BSC ccfDNA Plasma kit (Cat. #AS1480)	
Analyses:	qPCR	This protocol was developed by Promega Applications Scientists and is intended for research use only.
Sample Types:	Urine, plasma, cerebrospinal fluid (CSF)	Users are responsible for determining suitability of the protocol for their
Input:	1ml	application. Further information can be found in
Materials Required:	 Maxwell[®] RSC ccfDNA Plasma kit (Cat. #AS1480) 	Technical Manual #TM454, available at: <u>www.promega.com/protocols</u>

Maxwell[®] RSC Instrument (Cat. #AS4500)

or e-mail: techserv@promega.com

Protocol:

- 1. Sample preparation: Spin samples at 2,000 x g for 10 minutes at 4°C. Transfer cleared sample to a fresh tube and repeat the spin. NOTE: For best results do this step with fresh samples before freezing. A freeze/thaw step before centrifugation may result in cell lysis and release of gDNA into the cleared sample.
- 2. Add 1ml of cleared sample to Well #1 of the Maxwell[®] RSC ccfDNA plasma cartridge.
- 3. Add a plunger to Well #8 and 60μ l of Elution Buffer into each elution tubes.
- 4. Place cartridges and elution tubes in the Maxwell[®] RSC rack.
- 5. Run on Maxwell[®] RSC instrument with the RSC ccfDNA protocol.



Results: The above protocol was tested with 1ml of fresh urine, fresh plasma, and previous frozen CSF per DNA extraction (n=3).

Table 1: Concentrations of extracted ccfDNA determined by qPCR using 75bp amplicon primers.

	Mean (ng/µl)	STD
Plasma	0.064	0.008
Urine	0.027	0.001
CSF	1.897	0.177



Figure 1: qPCR ratios using 75bp and 300bp amplicon primers for analysis of quality of the ccfDNA. Ratios above 1 indicates smaller DNA fragments than longer fragments. Since ccfDNA is typically small, the larger the ratio indicates more ccfDNA and less gDNA.