

Bioluminescent Platform for Cellular Energy Metabolism Studies

Gediminas Vidugiris, Donna Leippe, Natasha Karassina, Mike Valley and Jolanta Vidugiriene


Promega Corporation, 2800 Woods Hollow Rd, Madison, WI 53711




1. Introduction

Energy uptake and utilization in eukaryotic cells is a dynamic process regulated by a series of interacting metabolic networks. Interrogation of this complex network relies on rapid, sensitive approaches that do not require extensive sample handling, are easily adaptable to 384-, 1536-well format and therefore can facilitate novel compound screening and evaluation. Utilizing a novel proluciferin substrate that in the presence of NAD(P)H is converted to luciferin and light production by luciferase, we developed a core technology for measuring key energy co-factors (NAD(P)/NAD(P)H) and extended it to multiple metabolite (glucose, lactate, glutamine, glutamate, glycerol, triglycerides) or enzyme activity assays. Here we demonstrate applicability of those assays to various formats (96-, 384- and 1536-well plates) and samples (2D, 3D cultures) using multiple automation platforms, including collecting the samples at different time points and performing multiple assays for the same sample using the Labcyte Echo acoustic dispenser.



2. Platform Overview: Based on NAD(P)H Detection




3. Applicability to Miniaturization and Automation




4. Multiple Metabolites from Small Amount of Cells with Minimal Sample Preparation



5. HTS Screening for Glycolysis and Glutaminolysis Inhibitors



6. Following Metabolite Consumption and Secretion During Cell Treatment



7. Conclusions

Bioluminescent energy metabolism platform is based on the core NAD(P)H detection technology:

- Provides sensitive detection (1-5 pmol/sample) with broad linear range (0.1-100µM) and wide dynamic range (max S/B >100 fold)
- Amendable to miniaturization with Z'>0.8 in 1536-well plates

The technology can be extended for multiple metabolite detection in different sample types:

- Glucose, Lactate, Glutamine, Glutamate, Glycerol, Triglyceride, Cholesterol, Cholesterol esters
- Cell medium and cells in culture, 3D models, tissues

Advantages:

- Multiple metabolites can be detected from the same sample
- Wide linear range provides flexibility when measuring metabolites over wide concentration range
- Optimized rapid in well cell lysis compatible with homogeneous format makes it suitable for HTS screening
- Multiplexing with viability/toxicity readouts provides normalization and allows separating effects on metabolism from global effects on cell health