

Cell-Based BRET Assays for Quantifying Antibody Mediated Blocking of Receptor-Ligand Interactions

Rachel Friedman Ohana¹, Robin Hurst¹, Mike Rosenblatt¹, Sergiy Levin², Cesear Corona², Thomas Machleidt¹, Matt Robers¹ and Keith V. Wood¹

¹Promega Corporation, 2800 Woods Hollow Rd, Madison, WI 53711; ²Promega Biosciences Inc, 277 Granada Dr, San Luis Obispo, CA 93401




1. Introduction

Antibody-based therapeutics have become a major trend in the drug discovery pipeline. Often these antibodies are generated to recognize particular ligands or receptors in order to block their interactions. To screen for such antibodies, we have developed a cell-based assay that rapidly quantifies ligand binding to a specific receptor on the cell surface. The assay utilizes bioluminescence resonance energy transfer (BRET) to detect the interaction of a fluorescently-labeled ligand with a receptor genetically tagged with a small and exceptionally bright NanoLuc® luciferase. Subsequently, the capacity of antibodies to block this interaction is quantified through the reduction in BRET signal.

Fluorescently-labeled ligands are commonly generated by random chemical modification of lysine residues. However, this approach routinely results in heterogeneous ligand populations that may exhibit variable interference in the interactions with their cognate receptors. We addressed this by developing an alternative labeling method that is both stoichiometric to minimize population variability and site specific to minimize perturbation of receptor binding. Using this method we generated three growth factors EGF, PDGF-B and VEGF-A₁₆₅ that were quantitatively labeled at the N-terminus and retained their biological function and binding affinity to their cognate receptors. In contrast, the random labeling resulted with significantly decreased bioactivity and binding affinity that was correlated to the number of labeled lysine residues.


Here, we used these site specifically labeled growth factors and their cognate receptor tyrosine kinases (RTK), which are important drug targets in cancer biology, to demonstrate the applicability of this cell-based assay for determining blocking efficacy of numerous therapeutic and research antibodies. This homogenous assay can be implemented as an early screening tool of biologics blocking efficacy and should be able to significantly advance antibody discovery work flows.

2. Measuring Antibody Blocking Efficacy by NanoBRET™




Highly specific: specificity governed by tight distance constraints (<5nM)
Signal strength: combination of NanoLuc® high emission intensity and a proper fluorophore
Live-cell assay: binding interactions are interrogated in a relevant cellular context
Homogeneous screen for blocking efficacy of therapeutic antibodies


3. Site Specific Labeling of Receptors Ligands




4. LC-MS/MS Analysis of Labeled Growth Factors




5. Influence of Labeling on Bioactivity




6. Influence of Labeling on Ligand Binding



7. Measuring Blocking Efficacy of Anti-Receptor Antibodies



8. Measuring Blocking Efficacy of Anti-Ligand Antibodies



9. Conclusions

Homogeneous NanoBRET™ assay to quantitate antibodies blocking efficacy

- **Cell-based assay:** cell surface receptor-ligand interactions are interrogated in a relevant cellular context
- **Rapid quantification:** direct measurements of antibodies capacity to block the interaction between a fluorescently labeled ligand and its cognate receptor that is genetically fused to NanoLuc®
- **NanoLuc®:** small and exceptionally bright
- **Labeled ligands:** labeled on a single N-terminal cysteine residue without any significant influence on their biological activity

The NanoBRET™ assay can be implemented as an early screening tool of biologics blocking efficacy