T7 RiboMAX™ Express RNAi System

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1. Elbashir, S.M. *et al.*(2001) RNA interference is mediated by 21- and 22-nucleotide RNAs. *Genes Dev.* **15**, 188-200.
2. Yang, D., Lu, H. and Erickson, J.W. (2000) Evidence that processed small dsRNAs may mediate sequence-specific mRNA degradation during RNAi in *Drosophila* embryos. *Curr. Biol.* **10**, 1191-200.
3. Fire, A. *et al.*(1998) Potent and specific genetic interference by double-stranded RNA in *Caenorhabditis elegans*. *Nature* **391**, 806-11.
4. Misquitta, L. and Paterson, B.M. (1999) Targeted disruption of gene function in *Drosophila* by RNA interference (RNA-i): A role for nautilus in embryonic somatic muscle formation. *Proc. Natl. Acad. Sci. USA* **96**, 1451-6.
5. Ngo, H. *et al.*(1998) Double-stranded RNA induces mRNA degradation in *Trypanosoma brucei*. *Proc. Natl. Acad. Sci. USA* **95**, 14687-92.
6. Sanchez-Alvarado, A. and Newmark, P.A. (1999) Double-stranded RNA specifically disrupts gene expression during planarian regeneration. *Proc. Natl. Acad. Sci. USA* **96**, 5049-54.
7. Lohmann, J.U. *et al.*(1999) Silencing of developmental genes in Hydra. *Dev. Biol.* **214**, 211-14.
8. Wargelius, A. *et al.*(1999) Double-stranded RNA induces specific developmental defects in zebrafish embryos. *Biochem. Biophys. Res. Commun.* **263**, 156-61.
9. Zamore, P.D. *et al.*(2000) RNAi: Double-stranded RNA directs the ATP-dependent cleavage of mRNA at 21 to 23 nucleotide intervals. *Cell* **101**, 25-33.
10. Hammone, S.M. *et al.* (2000) An RNA-directed nuclease mediates post-transcriptional gene silencing in *Drosophila* cells. *Nature* **404**, 293-6.
11. Schenborn, E.T. and Mierendorf, R.C. (1985) A novel transcription property of SP6 and T7 RNA polymerases: Dependence on template structure. *Nucl. Acids Res.* **13**, 6223-36.
12. Grishok, A., Tabara, H. and Mello, C.C. (2000) Genetic requirements for inheritance of RNAi in *C. elegans*. *Science* **287**, 2494-7.
13. Tabara, H., Grishok, A. and Mello, C.C. (1998) RNAi in *C. elegans*: Soaking in the genome sequence. *Science* **282**, 430-1.
14. Timmons, L. and Fire, A. (1998) Specific interference by ingested dsRNA. *Nature* **395**, 854.
15. Kamath, R.S. *et al.*(2001) Effectiveness of specific RNA-mediated interference through ingested double-stranded RNA in *Caenorhabditis elegans*. *Genome Biol.* **2**, RESEARCH0002.
16. Fraser, A.G. *et al.*(2000) Functional genomic analysis of *C. elegans* chromosome I by systematic RNA interference. *Nature* **408**, 325-30.
17. Clemens, J.C. *et al.*(2000) Use of double-stranded RNA interference in *Drosophila* cell lines to dissect signal transduction pathways. *Proc. Natl. Acad. Sci. USA* **97**, 6499-503.
18. Caplen, N.J. *et al.* (2000) dsRNA-mediated gene silencing in cultured *Drosophila* cells: A tissue culture model for the analysis of RNA interference. *Gene* **252**, 95-105.
19. Levashina, E.A. *et al.*(2001) Conserved role of a complement-like protein in phagocytosis revealed by dsRNA knockout in cultured cells of the mosquito, *Anopheles gambiae*. *Cell* **104**, 709-18.
20. Biyasheva, A. *et al.*(2001) Glue secretion in the *Drosophila* salivary gland: A model for steroid-regulated exocytosis. *Dev. Biol.* **231**, 234-51.
21. Krieg, P.A. (1990) Improved synthesis of full-length RNA probe at reduced incubation temperatures. *Nucl. Acids Res.* **18**, 6463.

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For Laboratory Use.

(a) U.S. Pat. No. 5,552,302, European Pat. No. 0 422 217, Australian Pat. No. 646803 and Japanese Pat. No. 3009458 have been issued to Promega Corporation for the methods and compositions for production of human recombinant placental ribonuclease inhibitor.

(b) U.S. Pat. Nos. 4,966,964, 5,019,556 and 5,266,687, Australian Pat. Nos. 616881 and 641261 and other pending and issued patents, which claim vectors encoding a portion of human placental ribonuclease inhibitor, are exclusively licensed to Promega Corporation.

(c) The method of recombinant expression of Coleopteraluciferase is covered by U.S. Pat. Nos. 5,583,024, 5,674,713 and 5,700,673.

(d) The RiboMAX™ Large Scale RNA Production System T7 (Cat.# P1300) is covered by U.S. Pat. No. 5,256,555 and is sold under a license from Ambion, Inc. It is intended for research use only. Parties wishing to use this product for other applications should contact Ambion, Inc.

(e) U.S. Pat. No. 4,766,072 has been issued to Promega Corporation for transcription vectors having two different bacteriophage RNA polymerase promoter sequences separated by a series of unique restriction sites into which foreign DNA can be inserted.

(f) U.S. Pat. Nos. 5,658,548 and 5,808,041 and Australian Pat. No. 689815 have been issued to Promega Corporation for nucleic acid purification on silica gel and glass mixtures. Other patents are pending.

(g) Licensed under one or more of U.S. Pat. Nos. 5,487,993 and 5,827,657 and European Pat. No. 0 550 693.

(h) U.S. Pat. No. 5,981,235 and Australian Pat. No. 729932 have been issued to Promega Corporation for methods for isolating nucleic acids using alkaline protease. Other patents are pending.

(i) Australian Pat. No. 730718 and Singapore Pat. No. 64532 have been issued to Promega Corporation for an improved filtration system and method. Other patents are pending.

(j) U.S. Pat. Nos. 5,283,179, 5,641,641, 5,650,289 and 5,814,471, Australian Pat. No. 649289, European Pat. No. 0 553 234 and Japanese Pat. No. 3171595 have been issued to Promega Corporation for a beetle luciferase assay method, which affords greater light output with improved kinetics as compared to the conventional assay. Other patents are pending.

(k) The PCR process is covered by patents issued and applicable in certain countries. Promega does not encourage or support the unauthorized or unlicensed use of the PCR process. Use of this product is recommended for persons that either have a license to perform PCR or are not required to obtain a license.

(l) U.S. Pat. No. 6,242,235 has been issued to Promega Corporation for polymerase stabilization by polyethoxylated amine surfactants. Other patents are pending.

(m) Patent Pending.

(n) For research use only. Not for human diagnostics or therapeutics or the commercial sale of transgenic animals. Use for any purpose other than research requires permission from Allele Biotechnology, Inc.
This product is covered under license from Carnegie Institution of Washington. Commercial use of this product may require a separate license from Carnegie.

© 2003 Promega Corporation. All Rights Reserved.

Flexi, pGEM, Riboprobe, RNasin and Wizard are trademarks of Promega Corporation and are registered with the U.S. Patent and Trademark Office. RiboMAX and siLenrGene are trademarks of Promega Corporation.

SYBR is a registered trademark of Molecular Probes, Inc.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.