

## Certificate of Analysis

### pFN31K *Nluc* CMV-neo Flexi® Vector:

Part No.	Size
N132A	20µg

**Description:** The pFN31K *Nluc* CMV-neo Flexi® Vector<sup>(a,b)</sup> is configured to facilitate simple, efficient transfer of the gene of interest into a vector designed for genetic attachment of NanoLuc® luciferase to the N-terminus of the protein of interest. This vector provides constitutive protein expression in mammalian cells using the human cytomegalovirus (CMV) immediate early enhancer/promoter. The vector can be used for both stable and transient gene expression. The vector backbone contains a neomycin phosphotransferase gene to allow selection in *E. coli* with kanamycin or in mammalian cell lines with neomycin.

The pFN31K *Nluc* CMV-neo Flexi® Vector contains the following features:

- A **CMV immediate early enhancer/promoter** for constitutive expression in mammalian cells.
- A **T7 RNA polymerase promoter** for cell-free expression of NanoLuc® fusion proteins.
- The **NanoLuc® gene** for high-sensitivity detection of NanoLuc® fusion proteins expressed in mammalian cells.
- The **lethal barnase gene** for positive selection of the insert. **Note:** The pFN31K *Nluc* CMV-neo Flexi® Vector can be propagated only in *E. coli* once the barnase gene is replaced with the protein-coding sequence of interest.
- A **neomycin/kanamycin-resistance gene** for selection in bacterial or mammalian cells.
- Unique **SgfI** and **PmeI** sites, which allow easy insertion of the sequence of interest from PCR product or compatible Flexi® vectors. When transferred in the proper reading frame, these sites create a readthrough sequence that can be joined to a protein-coding region flanked by SgfI and PmeI sites. In-frame transfer results in gene encoding a NanoLuc® fusion to the N-terminus of the protein of interest. Once inserted in this vector, the sequence is available for transfer to other Flexi® Vectors. For more information, see the *Flexi® Vector Systems Technical Manual* #TM254, available online at: [www.promega.com/protocols/](http://www.promega.com/protocols/)

**Concentration:** 100ng/µl.

**GenBank® Accession Number:** KF793053.

**Storage Buffer:** The pFN31K *Nluc* CMV-neo Flexi® Vector is supplied in 10mM Tris-HCl, 1mM EDTA (pH 8.0).

**Storage Conditions:** See Product Information Label for storage recommendations and expiration date. Avoid multiple freeze-thaw cycles and exposure to frequent temperature changes. These fluctuations can greatly alter product stability.

#### Usage Notes:

1. This vector was designed to be used with the Flexi® Vector System, a directional cloning method to shuttle protein-coding sequences between compatible vectors. To prepare the NanoLuc® fusion protein, the protein coding region is cloned into the pFN31K *Nluc* CMV-neo Flexi® Vector using the Flexi® System, Entry/Transfer (Cat.# C8640). For more information, see the *Flexi® Vector Systems Technical Manual* #TM254.

## Quality Control Assays

### Contaminant Assays

**Contaminating Nucleic Acids:** RNA, single-stranded DNA and chromosomal DNA are not evident in an overload sample of this vector as determined by agarose gel electrophoresis.

**Nuclease Assay:** Following incubation of 1µg of the vector in Restriction Enzyme Buffer at 37°C for 16–24 hours, no evidence of nuclease activity is detected by agarose gel electrophoresis.

**Physical Purity:**  $A_{260}/A_{280} \geq 1.80$ ,  $A_{260}/A_{250} \geq 1.05$ .

### Functional Assays

**Identity Assay:** The vector has been sequenced completely and has 100% identity with the published sequence available at: [www.promega.com/vectors/](http://www.promega.com/vectors/)

**Restriction Enzyme Digests:** Vector DNA is analyzed for the presence of certain restriction enzyme sites by incubation with a variety of restriction enzymes at the specified digestion temperature for 1 hour. Samples are examined by agarose gel electrophoresis, comparing cut and uncut vector DNA with marker DNA.

Part# 9PIN132

Printed 11/13



AF9PI N132 1113N132



**Promega**

#### Promega Corporation

2800 Woods Hollow Road	
Madison, WI 53711-5399	USA
Telephone	608-274-4330
Toll Free	800-356-9526
Fax	608-277-2516
Internet	<a href="http://www.promega.com">www.promega.com</a>

#### PRODUCT USE LIMITATIONS, WARRANTY, DISCLAIMER

Promega manufactures products for a number of intended uses. Please refer to the product label for the intended use statements for specific products. Promega products contain chemicals which may be harmful if misused. Due care should be exercised with all Promega products to prevent direct human contact.

Each Promega product is shipped with documentation stating specifications and other technical information. Promega products are warranted to meet or exceed the stated specifications. Promega's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Promega makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO PROMEGA PRODUCTS. In no event shall Promega be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of Promega products to perform in accordance with the stated specifications.

© 2013 Promega Corporation. All Rights Reserved.

Flexi, FuGENE, HaloTag, Nano-Glo and NanoLuc are registered trademarks of Promega Corporation. NanoBRET is a trademark of Promega Corporation.

GenBank is a registered trademark of U.S. Dept. of Health and Human Services.

All specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.

Part# 9PIN132  
Printed in USA 11/13.

Signed by:

*J. Stevens*

J. Stevens, Quality Assurance

## pFN31K *Nluc* CMV-neo Flexi® Vector Features

The following features are present based on nucleotide sequence.

CMV immediate early enhancer/promoter	1–742
Chimeric intron	857–989
T7 RNA polymerase promoter (–17 to +3)	1033–1052
NanoLuc® protein coding region	1065–1577
Linker	1578–1589
SgfI	1590–1597
Barnase coding region	1621–1956
PmeI	1958–1965
SV40 late poly(A) region	2117–2338
SV40 early enhancer/promoter	2437–2849
EM7 bacterial promoter	2863–2929
Neomycin phosphotransferase coding region	2943–3737
Synthetic poly(A) signal	3801–3849
ColE1-derived plasmid replication origin	4085–4121




Figure 1. pFN31K *Nluc* CMV-neo Flexi® Vector circle map.

## Related Products

Product	Size	Cat. #
Flexi® System, Entry/Transfer	5 entry and 20 transfer reactions	C8640
Flexi® System, Transfer	100 transfer reactions	C8820
Carboxy Flexi® System, Transfer	50 transfer reactions	C9320
10X Flexi® Enzyme Blend (SgfI & PmeI)	25µl	R1851
	100µl	R1852
Nano-Glo® Luciferase Assay	10ml	N1110
	100ml	N1120
	10 × 10ml	N1130
	10 × 100ml	N1150
FuGENE® HD Transfection Reagent	1ml	E2311
	5 × 1ml	E2312

©BY USE OF THIS PRODUCT, RESEARCHER AGREES TO BE BOUND BY THE TERMS OF THIS LIMITED USE LABEL LICENSE. If researcher is not willing to accept the conditions of this limited use statement, and the product is unused, Promega will accept return of the unused product and provide researcher with a full refund.

Researchers may use this product for research use only; no transfer or commercial use of this product is allowed. Commercial use means any and all uses of this product or derivatives by a party in exchange for consideration, including, but not limited to (1) use in further product manufacture; (2) use in provision of services, information or data; and (3) resale of the product or its derivatives, whether or not such product or derivatives are resold for use in research. Researchers shall have no right to modify or otherwise create variations of the nucleotide sequence of the luciferase gene except that researcher may (1) create fused gene sequences provided that the coding sequence of the resulting luciferase gene has no more than four deoxynucleotides missing at the affected terminus compared to the intact luciferase gene sequence; and (2) insert and remove nucleic acid sequences in splicing research predicated on the inactivation or reconstitution of the luminescence of the encoded luciferase. No other use of this product or derivatives is authorized without the prior express written consent of Promega.

### In addition, researchers must:

- (1a) use Nano-Glo®-branded luminescent assay reagents (LARs) for all determinations of luminescence activity of this product and its derivatives; or
- (1b) contact Promega to obtain a license for use of the product and its derivatives with LARs not manufactured by Promega.

**For uses of Nano-Glo®-branded LARs intended for energy transfer (such as bioluminescence resonance energy transfer) to acceptors other than a genetically encoded autofluorescent protein, researchers must:**

- (2a) use NanoBRET™-branded energy acceptors (e.g., BRET-optimized HaloTag® ligands) for all determinations of energy transfer activity by this product and its derivatives; or
- (2b) contact Promega to obtain a license for use of the product and its derivatives for energy transfer assays to energy acceptors not manufactured by Promega.

Researchers may transfer derivatives to others for research use provided that at the time of transfer a copy of this label license is given to the recipients and recipients agree to be bound by the terms of this label license. With respect to any uses outside this label license, including any diagnostic, therapeutic, prophylactic or commercial uses, please contact Promega for supply and licensing information. PROMEGA MAKES NO REPRESENTATIONS OR WARRANTIES OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING FOR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, WITH REGARD TO THE PRODUCT. The terms of this agreement shall be governed under the laws of the State of Wisconsin, USA.

©U.S. Pat. No. 8,557,970 and other patents pending.